


# Hard X-ray Diagnostics in the HSX

A. E. Abdou, A. F. Almagri, D. T. Anderson, J. N. Talmadge, S. P. Gerhardt, R. R. O'Connell\*

The HSX Plasma Laboratory, University of Wisconsin-Madison, USA

\* MST Plasma Laboratory, University of Wisconsin-Madison, USA

E-mail: [abdou@cae.wisc.edu](mailto:abdou@cae.wisc.edu)



APS, 44th Annual Meeting of the Division of Plasma Physics  
November 11-15, 2002; Orlando, Florida



THE UNIVERSITY  
of  
**WISCONSIN**  
MADISON



# Abstract

**A. E. Abdou, A. F. Almagri, D. T. Anderson, J. N. Talmadge, S. P. Gerhardt, R. R. O'Connell\***  
**The HSX Plasma Laboratory, University of Wisconsin-Madison, USA**  
**\* MST Plasma Laboratory, University of Wisconsin-Madison, USA**  
**E: [abdou@cae.wisc.edu](mailto:abdou@cae.wisc.edu)**

In the Helically Symmetric eXperiment (HSX), electrons are ECRH heated using 28 GHz, 2nd harmonic gyrotron. The normal configuration is Quasi-Helical Symmetric (QHS), it has a dominant  $n=4$ ,  $m=1$  components in the magnetic field spectrum. With a set of auxiliary coils, the quasihelical symmetry can be broken (MIRROR, HILL and WELL, modes). In this work the resolved hard X-ray emission in the HSX is analyzed using CdZnTe detector (1  $\mu$ s shaping time, Gain =50). The detector is housed in a lead box adjacent to the vacuum vessel, with 0.8 mm pinhole and 200  $\mu$ m SS filter placed in front of the detector. The hard X-ray pulse height spectrum was accumulated in a series of similar ECRH discharges. The behavior of the fast electrons has been studied for densities in the range of 0.2 to 0.9  $\times 10^{12}$  cm<sup>-3</sup>. The magnetic configuration has also been altered between QHS and MIRROR, and ANTI-MIRROR modes in order to determine the effect of magnetic ripple on characteristic energies and densities of fast electrons. Pulse height analysis of hard x-ray emission shows the presence of X-ray photons with energies as high as 650 KeV. Hard X-ray emission is strong function of plasma density; where at low density (0.2  $\times 10^{12}$  cm<sup>-3</sup>) the intensity is ~80 times higher than at higher density (0.9  $\times 10^{12}$  cm<sup>-3</sup>) for QHS mode. The hard X-ray intensity is higher in the QHS than in the MIRROR mode. Also the high energy tail of the spectrum extends to higher energies in QHS (~ 650 KeV). Future work will involve transport calculation to determine fast electron energy distribution.

\*This work is supported by US DOE grant number DE-FG02-93ER54222.



# Outline

- **Abstract**
- **HSX Device Parameters.**
- **Detector and Hardware Description.**
- **Source of fast electrons in HSX.**
- **Pulse Height Analysis program.**
- **HXR measurement results:**
  - Density Scan.**
  - Magnetic Configuration Scan.**
  - Neutral Pressure Scan.**
- **Conclusion.**




# Conclusion

- **Hard X-ray signal last for at least 50 ms after the ECRH is turned off indicating good confinement of fast electrons.**
- **Hard X-ray emission is a strong function of plasma density; where at low density ( $0.2 \times 10^{12} \text{ cm}^{-3}$ ) the intensity is  $\sim 80$  times higher than at high density ( $0.9 \times 10^{12} \text{ cm}^{-3}$ ).**
- **Hard Hard X-ray intensity is very low in ANTI-MIRROR mode at all plasma densities (poor confinement).**
- **The hard X-ray spectrum extends to energies as high as 650 keV in the QHS mode.**
- **Increasing the neutral base pressure beyond  $4.4 \times 10^{-6}$  suppress the fast electrons that are produced by toroidal voltage.**







# HSX Device Parameters

Major Radius	1.2 m
$\langle r \rangle$	0.11 m
Volume	$\sim .44 \text{ m}^3$
Field Periods	4
$i_{\text{axis}}$	1.05
$i_{\text{edge}}$	1.12
Coils/period	12
$B_0$ (max.)	1.25 T
Magnetic Field Flattop Length	0.2 s
Auxiliary Coils	48



# QHS, MIRROR and ANTI-MIRROR Modes

- QHS and Mirror have different magnetic field spectra...  
...but similar well depth, rotational transform, surface shape, and plasma volume.




APS, 44th Annual Meeting of the Division of Plasma Physics  
November 11-15, 2002; Orlando, Florida




# Sources of Fast Electrons in HSX

*HXR signal preionization off (Filament off)*



*HXR signal with preionization (Filament On)*




**Fast electrons in HSX are generated from two different mechanisms**

**1- dB/dt source**

**2- ECRH Heating source**



# Loop Voltage in HSX



- Loop Voltage time variation which gives rise to dB/dt source





# Detector

- **Detector Type:** CdZnTl
- **Detector Characteristics:**

**Dimension:** 10x10x2 mm

**Peak/Valley:**  $> 8:1$  @ 59.5 KeV

**Resolution:**  $< 10\%$  (6 KeV)@ 59.5 KeV (FWHM)

**Peak/Valley:**  $> 3:1$  @ 122 KeV

**Resolution:**  $< 6\%$  (8 KeV)@ 122 KeV (FWHM)



**Peak/Valley:**  $> 1.8:1$  @ 622 KeV

**Resolution:**  $< 3\%$  (20 KeV)@ 622 KeV (FWHM)




# Detector (continue)

- The detector has a built in pre-amplifier, it is connected to a high voltage power supply, and an Ortec EG&G 671 amplifier.
- The amplifier output signal is then connected to a data acquisition system.
- The raw signal is then processed using a pulse height analysis program.




# Detector Calibration

Number of Photons Vs Pulse Amplitude for Different Gain Settings Using Am241 (60 KeV)




Pulse Amplitude Vs Gain Using Am241 (60 KeV)



# Filter Transmission

X-ray Transmission for Fe from 0.001 to 2 Mev



$$I(x) = I_0 e^{-\mu x}$$

$$\text{Transmission} = I(x) / I_0 = e^{-\mu x}$$

Where:  $X = 200 \mu\text{m}$ ,  $\rho_{\text{Fe}} = 7.874 \text{ gm/cm}^2$  and  $\mu_{\text{Fe}}$  mass absorption coefficient for Fe ( $\text{cm}^2/\text{gm}$ )


Transmission is calculated from the above formula using tabulated values for  $\mu_{\text{Fe}}$  (from <http://physics.nist.gov/PhysRefData/FFast/html/form.html>)










# Pulse Height Analysis

The pulse height analysis program is written in IDL. The program mainly evaluate the following:

- 1- Resolving the Gaussian signals (single & double).
- 2- Least square fitting for the signals
- 3- Spectral Analysis of the signal
- 4- Calculation of the electron Temperature




# Pulse Heights Analysis (continue)




# HXR Spectrum at Low Density (QHS mode)

HXR Spectrum for QHS mode ( $n = 2 \times 10^{11} \text{ cm}^{-3}$ )




# HXR Spectrum at intermediate density (QHS mode)






# HXR Spectrum at High density (QHS mode)

HXR Spectrum for QHS mode ( $n = 9 \times 10^{11} \text{ cm}^{-3}$ )




# HXR Spectrum at Low density (MIRROR mode)

HXR Spectrum for Mirror mode ( $n = 2 \times 10^{11} \text{ cm}^{-3}$ )




# HXR Spectrum at Intermediate density (MIRROR mode)

HXR Spectrum for Mirror mode ( $n = 6 \times 10^{11} \text{ cm}^{-3}$ )




# HXR Spectrum at High density (MIRROR mode)




# HXR Spectrum at Low density (ANTI- MIRROR mode)

HXR Spectrum for Anti-Mirror mode ( $n = 2 \times 10^{11} \text{ cm}^{-3}$ )




# HXR Spectrum at Intermediate density (ANTI-MIRROR mode)

HXR Spectrum for Anti-Mirror mode ( $n = 5 \times 10^{11} \text{ cm}^{-3}$ )




# HXR Spectrum at High density (ANTI-MIRROR mode)

HXR Spectrum for Anti-Mirror mode ( $n = 9 \times 10^{11} \text{ cm}^{-3}$ )



# Summary

Number of Counts Vs Density for Different Magnetic Configuration




- Hard X-ray emission is strong function of plasma density; where at low density ( $0.2 \times 10^{12} \text{ cm}^{-3}$ ) the intensity is  $\sim 80$  times higher than at high density ( $0.9 \times 10^{12} \text{ cm}^{-3}$ ).
- The HXR spectrum extends to energies as high as 650 keV.






# Base Pressure Scan

*HXR signal (Base pressure =  $7.5 \text{ E-}08$ )*



*HXR signal (Base pressure =  $5.6 \text{ E-}07$ )*




**Increasing the tank base pressure decreases the the HXR intensity**




# Base Pressure Scan (Continue)

*HXR signal (Base pressure =  $6.5 \text{ E-07}$ )*



*HXR signal (Base pressure =  $4.4 \text{ E-06}$ )*



**Increasing the tank base pressure beyond  $4.4 \times 10^{-06}$  will suppress the fast electrons (No HXR signal)**



# Conclusion

- **Hard X-ray signal last for at least 50 ms after the ECRH is turned off indicating good confinement of fast electrons.**
- **Hard X-ray emission is a strong function of plasma density; where at low density ( $0.2 \times 10^{12} \text{ cm}^{-3}$ ) the intensity is  $\sim 80$  times higher than at high density ( $0.9 \times 10^{12} \text{ cm}^{-3}$ ).**
- **Hard Hard X-ray intensity is very low in ANTI-MIRROR mode at all plasma densities (poor confinement).**
- **The hard X-ray spectrum extends to energies as high as 650 keV in the QHS mode.**
- **Increasing the neutral base pressure beyond  $4.4 \times 10^{-6}$  suppress the fast electrons that are produced by toroidal voltage.**

